

# TASKI® Degragerm

Detergent, disinfectant and deodoriser for floors and surfaces

## Description

Highly concentrated daily cleaner for all types of water-resistant hard surfaces and floors, including protected floors.

## Key properties

- Bactericidal, fungicidal and virucidal activities
- Cleans without altering the shine of emulsions
- Fresh pine fragrance
- Low foam profile

## Benefits

- Disinfecting properties guarantees hygiene and safety by fighting against crosscontamination
- Can be used on most floors and surfaces, including protected floors
- Pine perfume increases the perception of cleanness in your premises, leaving a fresh fragrance
- Can be used in a scrubber drier as well as mop applications

## Use instructions

### Dosage:

Degragerm is available in various pack sizes and is either dosed manually (standard concentrate) or through one of our dosage equipments (J-Flex or QuattroSelect). For detailed information regarding the specific dosing equipment, please refer to the Technical data table. Proper dosage saves costs and minimizes environmental impact.

## Application

### Spray method:

- Spray the solution on a damp cloth and wipe
- Use a sponge pad to remove stubborn soil
- Rinse or replace the cloth regularly

### Bucket method:


- Apply the solution with a cloth/sponge or mop and wipe
- Use a sponge pad to remove stubborn soil

### Wet mopping method:

- Apply the solution with a mop and remove soiled solution

### Machine method:

- Fill in the solution into the full water tank of a single disc machine or scrubber drier using the high flow eductor through the filling gun.
- Apply the solution, allow if necessary and remove the soiled solution.



# TASKI<sup>®</sup> Degragerm

Detergent, disinfectant and deodoriser for floors and surfaces

## Technical data

Platform	SKU	Description/Pack size	Spray bottle	Bucket	Appearance	pH value neat	pH value in-use	Relative density (20°C)
Standard concentrate	7517609	TASKI Degragerm 6 x 1 L	5ml (1%) / 500 ml	50 ml (0,5%) / 10 L	Clear, green liquid	11	10.5	1.02 g/cm <sup>3</sup>
J-Flex	7517610	TASKI Degragerm J-Flex 1 x 1.5 L	300 x 500 ml (1%)	60 x 5 L (0,5%)	Clear, green liquid	11	10.5	1.05 g/cm <sup>3</sup>
SmartDose	7517843 7518810	TASKI Degragerm SD 1 x 1.4 L	350 x 500 ml (1%)	52 x 5 L (0,5%)	Clear, green liquid	11	10.5	1.05 g/cm <sup>3</sup>
QuattroSelect	7517568	TASKI Degragerm QS 2 x 2.5 L	1%	0,5 % (light soiling) 1 % (heavy soiling)	Clear, green liquid	11	10.5	1.05 g/cm <sup>3</sup>

The above data is typical of production and should not be taken as specification. For manual dosing of 1L, 5L products, corresponding dispenser pumps are available. Empty spray bottle and corresponding sprayers are also available.

## Safe handling and storage information

Full guidance on the handling and disposal of this product is provided in a separate Safety Data Sheet (SDS). <https://sds.sealedair.com/>. Store in original closed containers away from extremes of temperature.

## Product compatibility

Do not use on water-sensitive surfaces/materials (e.g. unsealed wood or cork). Test material compatibility on a small, inconspicuous place before use.

## Environmental information

The surfactants used in this product are biodegradable in compliance with the requirement of EU Detergent Regulation, EC 648/2004.

## Microbiological data

Bactericidal acc. EN1276 (0,5% / 5min, 20°C /dirty conditions);  
Yeasticidal acc. EN 1650 (0,5% / 5min, 20°C /dirty conditions);  
Fongicidal acc. EN1650 (0,5% / 5min, 20°C /dirty conditions);  
Active against Hepatitis B at (1% / 5 min) or (0.5% / 15 min / clean conditions);  
Active against Hepatitis B at (1.0% /15 min / dirty conditions);  
Active against HIV at (0.25% /5 min / clean conditions).